


# Painted Steel Casework



## Design Appearance:

- Our quick delivery cabinetry is offered in a flush inset style. All cabinet component parts of the unit are manufactured ensuring uniformity, interchangeability and accurate alignment.

## Color:

- Our quick delivery cabinetry is offered in a neutral color of Pearl White

## Field Convertible:

- Field convertible design allows for conversion of cabinet door and drawer front styles, drawer body suspension and cabinet horizontal and vertical rails, with the use of simple hand tools.

## LEED Compliance:

- Environmental compliant paints, glues and adhesives are specifically formulated to safeguard indoor air quality

## Industry Compliance:

- SEFA 8 Metal, NFPA 30, UL listed and FM approved

## Drawer Slide:

- 100 lbs. full extension, Standard

## Standard Features:

- Stainless Steel Wire Pull
- Stainless Steel Five Knuckle Hinge
- Nylon Roller Door Catch

## CornerStone Door and Drawer Size Code

Mod Code	Height
<b>Drawers</b>	
3M	3.649" (9.3 cm)
4M	4.909" (12.5 cm)
6M	7.429" (18.9 cm)
8M	9.949" (25.3 cm)
9M	11.209" (28.5 cm)
10M	12.469" (31.7 cm)
<b>Doors</b>	
15M	18.769" (47.7 cm)
18M	22.549" (57.3 cm)
20M	25.069" (63.7 cm)
22M	27.589" (70.1 cm)
24M	30.109" (76.5 cm)
28M	35.149" (89.3 cm)
37M	47.749" (121.3 cm)
63M	79.249" (201.3 cm)

## Painted Steel Casework

### **INSTOCK® Steel Laboratory Casework**

The **INSTOCK** program is ready when you are and is designed to ship complete laboratories within 10 days.


Our **INSTOCK** program provides you an extensive selection of lab furnishings that are available in a matter of days.\*

All **INSTOCK** steel cabinets and chemical hoods have a powder coated “Pearl White” finish that is SEFA 8 compliant.\*\* A wide range of other colors and cabinet styles are available in 6-8 weeks.



\*Most **INSTOCK**® items ship in 10 days. All items are subject to availability at the time of the order. Large quantity orders may require longer lead times.


\*\*For a color sample of Pearl White email: [orders@instocklabs.com](mailto:orders@instocklabs.com)



## CornerStone Sitting Height Base Cabinet


- Order work surfaces separately
- Provided without locks
- 27 9/16" (70 cm) Tall
- Standard Color: Pearl White
- All door and combination cabinets have a removable back
- All drawer cabinets have a cutout for service access

### Door Cabinets




JTP163-36	36" W (91.4 cm)
-----------	-----------------

### Drawer Cabinets




JTP174-18P	18" W (45.7 cm)
------------	-----------------




JTP070-24P	24" W (61 cm)
------------	---------------

### Door and Drawer Cabinets




Left Hinged	
JTP169-18P	18" W (45.7 cm)



## CornerStone ADA Height Base Cabinet


- Order work surfaces separately
- Provided without locks
- 32 5/8" (82.9 cm) Tall
- Standard Color: Pearl White
- All door and combination cabinets have a removable back
- All drawer cabinets have a cutout for service access

### Drawer Cabinets



JTP925-18P	18" W (45.7 cm)
------------	-----------------


### Sink Cabinet



Assembly Required

JTP981ACC-36	36" W (91.4 cm)
--------------	-----------------


# Painted Steel Casework



## CornerStone Standing Height Base Cabinet


- Order work surfaces separately
- Provided without locks
- 35 1/8" (89.2 cm) Tall
- Standard Color: Pearl White
- All door and combination cabinets have a removable back
- All drawer cabinets have a cutout for service access

### Door Cabinets




JTP103-30	30" W (76.2 cm)
JTP103-36	36" W (91.4 cm)
JTP103-48	48" W (121.9 cm)

### Drawer Cabinets




JTP122-18P	18" W (45.7 cm)
JTP122-24P	24" W (61 cm)




JTP121-36P	36" W (91.4 cm)
------------	-----------------

### Drawer Cabinets




JTP124-48P	48" W (121.9 cm)
------------	------------------

### Door/Drawer Cabinets




JTP106-36P	36" W (91.4 cm)
JTP106-42P	42" W (106.7 cm)
JTP106-48P	48" W (121.9 cm)




Right Hinged	
JTP118-18P	18" W (45.7 cm)
JTP118-24P	24" W (61 cm)
Left Hinged	
JTP119-18P	18" W (45.7 cm)
JTP119-24P	24" W (61 cm)

### Corner/Sink Cabinet




Right Hinged	
JTP190-32SE	32-7/8" W (83.5 cm)



## CornerStone Standing Height Base Cabinet


- Order work surfaces separately
- Provided without locks
- 35 1/8" (89.2 cm) Tall
- Standard Color: Pearl White
- All door and combination cabinets have a removable back
- All drawer cabinets have a cutout for service access

## Sink Cabinets



JTP181-36	36" W (91.4 cm)
JTP181-48	48" W (121.9 cm)

## End Island Sink Cabinets




JTP186-52	52" W (132.1 cm)
JTP186-58	58" W (147.3 cm)

## Sink Supports



JP535-36	36" W (91.4 cm)
JP535-48	48" W (121.9 cm)
JP535-52	52" W (132.1 cm)
JP535-58	58" W (147.3 cm)


## Acid Storage Cabinets



- Full-Height doors are lined with 1/8" thick polyethylene
- Chemical resistant 3/16" thick molded Polyethylene lining with coved corners and 2" lip at front of cabinet opening
- Includes vent kit
- One half-depth shelf

JTP185-24GG (Left Hinged)	24" W (61 cm)
JTP188-36GG	36" W (91.4 cm)
JTP188-48GG	48" W (121.9 cm)


## Solvent Storage Cabinets



- Construction meets requirements listed by UFC, OSHA, and NFPA No. 30-1993
- Two 2" pipe threaded vent outlets on back
- Rear and top access panels
- Lock is included
- Self-closing


JTP189-24SLAS (Right Hinged)	24" W (61 cm)
JTP189-36SLASY	36" W (91.4 cm)
JTP189-48SLASY	48" W (121.9 cm)

## Standard Hood Cabinets



JTP103-30	30" W (76.2 cm)
JTP103-36	36" W (91.4 cm)
JTP103-48	48" W (121.9 cm)


# Painted Steel Casework



## CornerStone Mobile Cabinet


- Order work surfaces separately
- Provided without locks
- 27 1/4" Tall and 22" Deep
- Standard Color: Pearl White
- Cabinets include 4 each heavy duty 3" casters, the front 2 casters are locking
- Drawer interlock system included and counter weighted anti-tip design.

### Combination Cabinets




Right Hinged	
JTP1068-18	18" W (45.7 cm)
JTP1068-24	24" W (61 cm)

### Drawer Cabinets




JTP1072-18	18" W (45.7 cm)
JTP1072-24	24" W (61 cm)



## CornerStone Mobile Cabinet


- Order work surfaces separately
- Provided without locks
- 32 3/8" Tall and 22" Deep
- Standard Color: Pearl White
- Cabinets include 4 each heavy duty 3" casters, the front 2 casters are locking
- Drawer interlock system included and counter weighted anti-tip design.

### Combination Cabinets




Right Hinged	
JTP1018-18	18" W (45.7 cm)
JTP1018-24	24" W (61 cm)

### Drawer Cabinets




JTP1021-18	18" W (45.7 cm)
JTP1021-24	24" W (61 cm)



## Mobile Cabinet Work Surface

- Work surface is 3/4" thick black phenolic resin that provide a 1" overhang on the front and a 1/2" overhang at the rear and sides.
- It adds 3/4" to the overall height of the cabinet.
- There is a grip notch included under the front edge of the work surface to help maneuver cabinet


Part Number	For use with:
MWS-18	18" (45.7 cm) wide mobile cabinet
MWS-24	24" (61 cm) wide mobile cabinet




## CornerStone Suspended Cabinet

- 23" Tall and 22" Deep
- Standard Color: Pearl White
- Designed to be hung from any of our table systems
- Each suspended cabinet is furnished with complete hanging hardware.


## Drawer Cabinets



JTP820-18	18" W (45.7 cm)
JTP820-24	24" W (61 cm)



JTP822-15	15" W (38.1 cm)
-----------	-----------------




PPD1822	18" W (45.7 cm)
---------	-----------------

## CornerStone Suspended Pencil Drawer

- 3 7/8" Tall and 22" Deep
- Standard Color: Pearl White
- Designed to be hung from any of our table systems
- Each suspended cabinet is furnished with complete hanging hardware.


# Painted Steel Casework



## CornerStone Apron Frame with Drawer

- Apron frame is 4" (11.4 cm) tall
- Use to span knee space area between two cabinets
- Attachment hardware included
- Standard color: Pearl White

Part Number	Width	Depth
JTPTBA-3022-1P	30" (76.2 cm)	22" (55.9 cm)
JTPTBA-3622-1P	36" (91.4 cm)	22" (55.9 cm)




## CornerStone Apron Rail

- Apron rails are 2" tall and 2" deep
- Used to span and support knee space areas between two cabinets
- Attachment hardware is included
- Standard color: Pearl White
- Can be cut to fit in the field



Part Number	Width
JPTBR-3002	30" (76.2 cm)
JPTBR-3602	36" (91.4 cm)
JPTBR-4802	48" (121.9 cm)
JPTBR-7202	72" (182.9 cm)

## Example use of Rails




## Panel End Supports

- Use in place of a cabinet at one end of a cabinet run to support drawer frames or apron rails
- Can be used to support the intersection of two aprons or rails
- Each panel has 2 leveling feet and is 2" (5.1 cm) wide
- Standard Color: Pearl White



Part Number	Depth	Height
<b>Sitting Height</b>		
JPTPO-2822	22" (55.9 cm)	27 9/16" (70 cm)
<b>ADA Height</b>		
JPTPO-3222	22" (55.9 cm)	32 5/8" (82.8 cm)
<b>Standing Height</b>		
JPTPO-3522	22" (55.9 cm)	35 1/8" (89.2 cm)




## Heavy Duty Wall Mount Bracket Countertop Support

- Used to support spans of work surface without structure mounted on the floor
- Adequate in-wall blocking is required.
- 11 gauge bracket has a 2" wide return flanges with mounting holes at the top and rear of the bracket
- Bracket is 24" tall and can be mounted at any chosen height
- Standard color: Pearl White

Part Number	Depth
FF-CTOP-18	18" (45.7 cm)
FF-CTOP-23	23" (58.4 cm)


# Painted Steel Casework



## CornerStone Wall Cabinets


- Provided without locks
- 30 5/8" (77.8 cm) Tall x 13" (33 cm) Deep
- Standard Color: Pearl White
- Two adjustable shelves
- Flush Bottom

## Solid Door Wall Cabinets




JTP203-36	36" W (91.4 cm)
JTP203-48	48" W (121.9 cm)

## Framed Glass Wall Cabinets



JTP206-36	36" W (91.4 cm)
JTP206-48	48" W (121.9 cm)

## Rear Wall Cabinet Filler




JC1392N0124P
--------------

## CornerStone Steel Solid Door Tall Cabinet

- Five shelves are included and are adjustable in half-inch increments
- Must be anchored to wall
- 84-5/16" Tall x 22" Deep
- Standard Color: Pearl White
- Tall cabinet with no lock shown below

### Solid Door Tall Cabinets




Tall Cabinet with No Lock	
JTP303-36	36" W (91.4 cm)
JTP303-48	48" W (121.9 cm)
Tall Cabinet with Standard Lock	
JTP303-36L	36" W (91.4 cm)
JTP303-48L	48" W (121.9 cm)
Tall Cabinet with Three-Point Latch	
JTP303-363PL	36" W (91.4 cm)
JTP303-483PL	48" W (121.9 cm)

## CornerStone Steel Glass Door Tall Cabinet

- Five shelves are included and are adjustable in half-inch increments
- Must be anchored to wall
- 84-5/16" Tall x 22" Deep
- Standard Color: Pearl White
- Tall cabinet with standard lock shown below


### Framed Glass Door Tall Cabinets



Tall Cabinet with Standard Lock	
JTP307-36L	36" W (91.4 cm)
JTP307-48L	48" W (121.9 cm)
Tall Cabinet with Three-Point Latch	
JTP307-363PL	36" W (91.4 cm)
JTP307-483PL	48" W (121.9 cm)




# Painted Steel Casework

## Knee Space Panels for use with Drawer Aprons



Part Number	Width
<b>Standing Height</b>	
JKSH-3004GG	30" (76.2 cm)
JKSH-3604GG	36" (91.4 cm)

## Knee Space Panels for use with 2" High Apron Rails




Part Number	Width
<b>Sitting Height</b>	
JKSL-ADJGG*	ADJUSTABLE 12"-48" (30.5 - 121.9 cm)

Part Number	Width
<b>ADA Height</b>	
JKSA-ADJGG*	ADJUSTABLE 12"-48" (30.5 - 121.9 cm)

Part Number	Width
<b>Standing Height</b>	
JKSH-30GG	30" (76.2 cm)
JKSH-36GG	36" (91.4 cm)
JKSH-ADJGG*	ADJUSTABLE 12"-48" (30.5 - 121.9 cm)
JKSH-72ADJGG*	ADJUSTABLE 12"-72" (30.5 - 182.9 cm)

\* Cut to fit and includes slip fit trim angle


## Corner Fillers



Part Number	Width
<b>Standing Height - 35 1/8"</b>	
JFCH-01	1" (2.5 cm)
JFCH-ADJ*	ADJUSTABLE 2"-12" (5.1 - 30.5 cm)

\* Cut to fit and includes slip fit trim angle


## Front Fillers Panels



- For use between adjacent base cabinets or between a cabinet and the wall

Part Number	Width
<b>Sitting Height - 27 9/16"</b>	
JFFL-ADJ*	2"-18" (5.1 - 45.7 cm)
<b>Standing Height - 35 1/8"</b>	
JFFH-01	1" (2.5 cm)
JFFH-04	4" (10.2 cm)
JFFH-08	8" (20.3 cm)
JFFH-ADJ*	2"-18" (5.1 - 45.7 cm)


## Rear and Center Fillers Panels



- For use between back-to-back base cabinets or to close space at the rear of the cabinet

Part Number	Width
<b>Sitting Height - 27 9/16"</b>	
JFRL-07	7" (18.7 cm)
<b>ADA Height - 32 5/8"</b>	
JFRA-07	7" (18.7 cm)
<b>Standing Height - 35 1/8"</b>	
JFRH-01	1" (3.5 cm)
JFRH-07	7" (18.7 cm)
JFRH-08	8" (21.3 cm)
JFRH-14	14" (36.5 cm)
JFRH-ADJ*	2" - 18" (5.1 - 45.7 cm)

## Pipe Enclosures - Steel



- Used to provide a cover for pipes running from the ceiling to the work surface
- Can be used for island or perimeter locations
- Can support pegboards or water purification systems
- Cover is pre-punched for attachment and is furnished in three sections for ease of access after installation

Part Number	Dimensions (W x D x H)
JP509-1284	12" x 6" x 84" (30.5 x 15.2 x 212.1 cm)
JP509-1884	18" x 6" x 84" (45.7 x 15.2 x 212.1 cm)